



PatientScript — Handwriting Capture

Profile: daniel-test-011


Page 1 of 2

1. Use a black or dark blue pen — gel or fine ballpoint works best.
2. Each letter gets 3 cells; punctuation gets 1. Let your handwriting vary naturally.
3. The horizontal lines between cells are your baseline. Sit each letter on it; descenders (g, j, p, q, y) dip below.
4. When done, photograph each page on a flat surface in good light. All four corner markers must be visible.

HOW THE BASELINE WORKS

The line runs BETWEEN cells, not through them.

Letters sit on the imaginary continuation. g j p q y dip below.



A	B	C	D	E					
F	G	H	I	J					
K	L	M	N	O					
P	Q	R	S	T					
U	V	W	X	Y					
Z									

PUNCTUATION

One sample of each. Periods and commas sit near the bottom; ' floats high.

.	,	'	-	!	?	:	;
---	---	---	---	---	---	---	---





Profile: daniel-test-011

Page 2 of 2

Sign your name once below. This appears at the end of each postcard.